

Fmoc-L-Lys(Dde)-OH

Cat. No.: HY-W041987 CAS No.: 150629-67-7 Molecular Formula: $C_{31}H_{36}N_2O_6$ Molecular Weight: 532.63

Amino Acid Derivatives Target:

Pathway: Others

Storage: Powder -20°C 3 years

2 years

In solvent -80°C 6 months

> -20°C 1 month

Product Data Sheet

SOLVENT & SOLUBILITY

In Vitro

DMSO: 100 mg/mL (187.75 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	1.8775 mL	9.3874 mL	18.7748 mL
	5 mM	0.3755 mL	1.8775 mL	3.7550 mL
	10 mM	0.1877 mL	0.9387 mL	1.8775 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (4.69 mM); Clear solution

BIOLOGICAL ACTIVITY

Description	Fmoc-L-Lys(Dde)-OH is a lysine derivative ^[1] .
In Vitro	Amino acids and amino acid derivatives have been commercially used as ergogenic supplements. They influence the secretion of anabolic hormones, supply of fuel during exercise, mental performance during stress related tasks and prevent exercise induced muscle damage. They are recognized to be beneficial as ergogenic dietary substances ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-1144.

 $\label{lem:caution:Product} \textbf{Caution: Product has not been fully validated for medical applications. For research use only.}$

Tel: 609-228-6898 Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Page 2 of 2 www.MedChemExpress.com